

VWR[®] FOR YOUR SPACE

Laboratory chairs and stools

Specialist seating

Office and visitor chairs

Other handy stuff

Are you sitting
comfortably?



VWR® FOR YOUR SPACE

VWR ensures you're sitting comfortably!

Choice – colour, design and type - without being overwhelming

Performance – market leading consistency and reliability backed up with warranties that run to many years

Value for money – well built products that last and don't cost a fortune to buy



INDEX

Seating

- General lab chairs and stools3
- Specialist chairs for cleanrooms and ESD.....11
- Office, other chairs and footrests.....13

Workspace products

- PCR Workstation15
- Granite anti-vibration table and slab15
- Lights16

Other handy stuff

- Steps and ladders17
- Trolleys18
- Boxes, organisers and holders.....23
- Log books and pens27
- Wipes and other cleaning materials incl. silicone grease29
- Ergonomic flooring33

THESE SYMBOLS INDICATE IMPORTANT PRODUCT FEATURES:

SYMBOLS



Warranty



Warranty



ESD



Warranty



Warranty

Chairs, Stools, Economy Line



132-0410

132-0411

132-0412

132-0413

132-0414

132-0415

132-0499

132-0500

All chairs and stools are available with a star base which is made from high strength sectional steel tubing, and is particularly robust and durable. The flat design prevents trip hazards.

- Anatomically-designed seat and backrest
- Washable, easy to clean, robust and durable
- Seat height adjustment by gas-lift

Glides: The glides are wear-proof plastic and have a large surface area and offering a high level of stability.

Castors: Dual-wheel castors with load-sensitive brakes have to be right for the floor. All chairs are supplied with soft castors for hard floors as standard. Hard castors for soft floors are also available as an option.

Stop&Go technology: For safety reasons chairs whose seat is more than 650 mm high must not be movable whilst in use. With innovative Stop&Go castors these chairs are only mobile when no one is sitting on them.

As well as the need for maximum hygiene and easy cleaning, laboratory chairs also have to meet a number of other requirements associated with routine laboratory tasks: Flexibility, in terms of the work applications and they must not take too much space.

Description	Height (mm)	Foot form	Foot ring	Pk	Cat. No.
Laminate wood chairs					
Laboratory chair, made from laminated beech, backrest adjustable in height and inclination, backrest pendular embedded	440 - 610	Glides	-	1	132-0410
Laboratory chair, made from laminated beech, backrest adjustable in height and inclination, backrest pendular embedded	440 - 610	Castors	-	1	132-0414
Laboratory chair, made from laminated beech, backrest adjustable in height and inclination, backrest pendular embedded	570 - 850	Glides	+	1	132-0412
Laboratory chair made from laminated beech, with the backrest adjustable in height and inclination	610 - 910	Stop&Go castors	+	1	132-0499
Artificial leather chairs					
Laboratory chair, made from artificial black leather, backrest adjustable in height and inclination, backrest pendular embedded	440 - 610	Glides	-	1	132-0411
Laboratory chair, made from artificial black leather, backrest adjustable in height and inclination, backrest pendular embedded	440 - 610	Castors	-	1	132-0415
Laboratory chair, made from artificial black leather, backrest adjustable in height and inclination, backrest pendular embedded	570 - 850	Glides	+	1	132-0413
Laboratory chair, made from artificial black leather, with the backrest adjustable in height and inclination	610 - 910	Stop&Go castors	+	1	132-0500

Description	Pk	Cat. No.
Accessories		
Castors	1 SET	132-0427

Laboratory chairs and stool, blue line



132-1375



132-1376



132-1379

These robust chairs and stool are designed for use in laboratories and are perfect for many applications from pharmaceutical, biotech, chemicals, health care and clean room to medical technology. The seat is made from artificial and antibacterial leather and the backrest is made of plastic in such a way that there are no joints. This means that the unit is extremely easy to keep clean.

- 3-D flex function on the front edge of the seat and the backrest follow automatically the movement of the body
- Solid plastic star base in black
- Coloured closed seat shell
- Easy to clean, washable and resistant to disinfectants
- Joint-free thanks to hygienic design
- Integrated handle ensures maximum mobility and simple handling
- Gas lift for quick height adjustment in sitting position
- Protective cover protects the gas lift from dust and dirt

Suitable for cleanrooms category 3 in accordance with ISO 14644-1.

Suitable for use in biotechnology laboratories: safety classes S1, S2, and S3.

GMP conformity. GS-certified.


The identification clip enables the chair to be easily, individually and quickly labelled. This means that every chair can easily be allocated to a user or an area. This create structure in the workplace. It can also be added and removed any time. Identification clip is included in a standard delivery.

Description	Height (mm)	Foot form	Foot ring	Pk	Cat. No.
Standard laboratory chairs					
Laboratory chair with seat's upholstery made blue artificial leather, plastic backrest and anthracite coloured closed seat shell	450–650	Soft castors for hard floors	-	1	132-1376
High laboratory chairs					
Laboratory chair with seat's upholstery made from blue artificial leather, plastic backrest and blue coloured closed seat shell	560–810	Stop&Go castors	Chrome plated, continiously height adjustable	1	132-1377
Laboratory chair with seat's upholstery made blue artificial leather, plastic backrest and anthracite coloured closed seat shell	560–810	Stop&Go castors	Chrome plated, continiously height adjustable	1	132-1378
Laboratory stool					
Laboratory stool with large seating surface, the seat upholstery is made from black PU-foam	450–650	Soft castors for hard floors	-	1	132-1379
Standard laboratory chairs					
Laboratory chair with seat's upholstery made from blue artificial leather, plastic backrest and blue coloured closed seat shell	450–650	Soft castors for hard floors	-	1	132-1375

Equipment Management powered by VEM technology


A modern Equipment Management system for the entire organisation

VEM - A complete solution to manage and monitor equipment in the lab



Cloud-based services

VEM delivers to the desktop. It is built using the latest and most powerful internet technologies to better support your requirements. VEM offers a modern widget-based user-interface that gives users access to all major functionality immediately.



- Centralised equipment repository
- Tracks all equipment service management needs: Service records, preventive maintenance schedules, creates service requests, etc.
- Enables strategic procurement decisions based on information stored: Total cost of ownership
- Contract and warranty management
- Runs reports and creates key performance indicators
- Compares performance of manufacturers, service providers, vendors, etc.

VEM combines your equipment portfolio, maintenance scheduling, inventory and work order management together for quick and easy administration, tracking and decision making.



For more information, visit our website vwr.com/VEM
Or send an email to: vwr catalyst.eu@vwr.com



132-0405

Chairs, stool, All-In-One line



All chairs and stools are available with a star base which is made from high strength sectional steel tubing, which is particularly robust and durable. The flat design prevents trip hazards.

- Complete hygienic design, easy to clean and disinfect
- Seat height adjustment by gas-lift
- Seat tilt adjustment
- Backrest adjustable in height and inclination

Glides: The glides are wear-proof plastic and have a large surface area offering a high level of stability.

Castors: Dual-wheel castors with load-sensitive brakes have to be right for the floor. All chairs are supplied with soft castors for hard floors as standard. Hard castors for soft floors are also available as an option.

Stop&Go technology: For safety reasons chairs whose seat is more than 650 mm high must not be movable whilst in use. With innovative Stop&Go castors these chairs are only mobile when no one is sitting on them.

SoftTouch PU foam, resistant against aggressive liquids. The permanent contact backrest follows all movements automatically and supports dynamic sitting, ergonomical shaped backrest for ideal fitting to the back.

Conform to security level S2, GMP

Description	Height (mm)	Foot form	Foot ring	Pk	Cat. No.
Laboratory chair, with SoftTouch PU foam, black	430 - 600	Glides	-	1	132-0405
Laboratory chair, with SoftTouch PU foam, black	430 - 600	Castors	-	1	132-0406
Laboratory chair, with SoftTouch PU foam, black	580 - 850	Glides	+	1	132-0407
Laboratory chair SoftTouch PU foam	600 - 870	Stop&Go castors	+	1	132-0498

Description	Pk	Cat. No.
Accessories		
Armrest, black	1 PAIR	132-0423
Castors	1 SET	132-0427
Height-adjustable foot ring for All-In-One chairs	1	139-0215

Chairs, All-In-One Alu



The VWR All-In-One Alu is the upgrade version of the All-In-One. All the chairs are designed with a polished aluminium 5 star base. The seat height adjustment allows the individual adjustment of the seat height depending on the user's seating habits, activity and height. Also, the seat inclination adjustment enables inclination of the entire seat up to 5°, providing supported active sitting – ensuring a vertical position of the pelvis, helping the spine to maintain its natural S shape.

- SoftTouch PU foam, resistant against aggressive liquids, easy to clean and disinfect
- Permanent contact backrest responds to any movement and supports dynamic sitting
- Seat height adjustment by gas-lift
- Ergonomically shaped backrest for ideal fitting of the back

Glides: The glides are wear-proof plastic and have a large surface area offering a high level of stability.

Castors: Dual-wheel castors with load-sensitive brakes have to be right for the floor. All chairs are supplied with soft castors for hard floors as standard. Hard castors for soft floors are also available as an option.

Stop&Go technology: For safety reasons chairs whose seat is more than 650 mm high must not be movable whilst in use. With innovative Stop&Go castors these chairs are only mobile when no one is sitting on them.

Conform to security level S2, GMP.



132-0746

Description	Height (mm)	Foot form	Foot ring	Pk	Cat. No.
All-In-One Alu polished aluminium chair	430 - 600	Glides	-	1	132-0743
All-In-One Alu polished aluminium chair	430 - 600	Castors	-	1	132-0744
All-In-One Alu polished aluminium chair	580 - 850	Glides	+	1	132-0745
All-In-One Alu polished aluminium chair	600 - 870	Stop&Go castors	+	1	132-0746

Chairs, stools, Mr. Lab



132-0416



132-0417

All chairs and stools are available with a five star base which is made from aluminium.

- Chair backrest adjustable in height and inclination
- Seat height adjustable by gas-lift
- Aluminium 5 star base

The perfect combination of comfort and ergonomics for optimal performance. Anatomically-moulded cushions with pelvis support and integrated lumbar support. Very comfortable, large and soft upholstery in artificial leather covering, washable and disinfectant-proof.

Glides: The glides are wear-proof plastic and have a large surface area offering a high level of stability.

Castors: Dual-wheel castors with load-sensitive brakes have to be right for the floor. All chairs are supplied with soft castors for hard floors as standard. Hard castors for soft floors are also available as an option.

Stop&Go technology: For safety reasons chairs whose seat is more than 650 mm high must not be movable whilst in use. With innovative Stop&Go castors these chairs are only mobile when no one is sitting on them.

Conform to security level S2, GMP

Description	Height (mm)	Foot form	Foot ring	Pk	Cat. No.
Chairs					
Laboratory chair, made from blue skai, with permanent contact backrest and aluminium polished base, plastic parts in black	470 - 620	Castors	-	1	132-0416
Laboratory chair, made from blue skai, with permanent contact backrest and aluminium polished base, plastic parts in black	470 - 620	Stop&Go castors	-	1	132-0939
Laboratory chair, made from blue skai, with permanent contact backrest and aluminium polished base, plastic parts in black	600 - 870	Glides	+	1	132-0417
Laboratory chair, made from blue skai, with permanent contact backrest and aluminium polished base, plastic parts in black	640 - 910	Stop&Go castors	+	1	132-0501

Description	Height (mm)	Foot form	Foot ring	Pk	Cat. No.
Chairs					
Laboratory chair, made from black skai, with permanent contact backrest and aluminium polished base, plastic parts in black	470 - 620	Soft castors	-	1	132-0965
Laboratory chair, made from black skai, with permanent contact backrest and aluminium polished base, plastic parts in black	600 - 870	Large scaled glides	Height adjustable, chrome plated	1	132-0966
Laboratory chair, made from black skai, with permanent contact backrest and aluminium polished base, plastic parts in black	620 - 890	Stop&Go castors	Height adjustable, chrome plated	1	132-0967
Stools					
Laboratory stool, plastic black frame, seat made from black skai; height adjustment based on pneumatic spring (plastic and steel part are black)	460 - 630	Castors	-	1	132-0425
Laboratory stool, plastic black frame, seat made from black PU foam; height adjustment based on pneumatic spring (plastic and steel part are black)	460 - 630	Castors	-	1	132-0426
Laboratory stool, extra large seat (Ø 400 mm) made from black skai leather, PA base (plastic and steel part are black)	600 - 870	Stop&Go castors	Height adjustable, chrome plated	1	132-0963
Laboratory stool, extra large seat (Ø 350 mm) made from PU foam, PA base (plastic and steel part are black)	600 - 870	Stop&Go castors	Height adjustable, chrome plated	1	132-0964
Description				Pk	Cat. No.
Accessories					
Armrest, black				1 PAIR	132-0423
Castors				1 SET	132-0427
Height-adjustable footring for 132-0416				1	139-0214



Chairs, Flex



Flex is the innovative solution for sitting and standing workplaces. Thus, Flex is usable both as stool and as standing rest, it is light, uncomplicated and easy to use. Flex adapts to its user and the workplace situation.


- Ergonomic, soft and comfortable seat with back support
- Turns through 360° for maximum freedom of movement
- Extra large seat height adjustment range by gas-lift
- SoftTouch PU foam, easy to clean and disinfect
- Polished aluminium 5 star base

Glides: The glides are wear-proof plastic and have a large surface area offering a high level of stability.

Castors: Dual-wheel castors with load-sensitive brakes have to be right for the floor. All chairs are supplied with soft castors for hard floors as standard. Hard castors for soft floors are also available as an option.

Stop&Go technology: For safety reasons chairs whose seat is more than 650 mm high must not be movable whilst in use. With innovative Stop&Go castors these chairs are only mobile when no one is sitting on them.

Description	Height (mm)	Foot form	Foot ring	Pk	Cat. No.
Flex lab standing rest	450 - 650	Castors	-	1	132-0747
Flex lab standing rest	510 - 780	Glides	-	1	132-0748
Flex lab standing rest	510 - 780	Glides	+	1	132-0749
Flex lab standing rest	540 - 810	Stop&Go castors	+	1	132-0750



VWR CATALYST
We Enable Science Through Services

Calibration and certification

Servicing and repairs

Maintenance contracts

Please visit us at: vwr.com/vwrcatalyst

Ergonomic laboratory chairs and stools



VWR132-1105



132-1099



132-1301

Nobody has a perfectly symmetric morphology. This is why a majority of people have back problems, sometimes with a feeling of heavy legs. The answer to this problem lies in the main features of this seat. Above all, a good working position is a dynamic posture that facilitates the bloodstream. The features of the seat help to re-balance the body naturally and improve blood flow in the legs. VWR® Ergonomic lab chairs and stools adapt automatically to the body and the movements of the user.

- Automatically adapts to the user
- Facilitates an upright posture
- Minimises spinal stress
- Seat height adjustment by gas lift
- Black steel, 5-branch, scuff resistant reinforced base
- Resistant against sparks, chemicals and weak acids

The innovative seat design consist of a short seat tilted forward, split into 2 parts, based on blue silent-blocks.

The PU foam seat has a thickness of 40 mm and great deformation capacity which provides better distribution of pressure when moving.

Description	Height (mm)	Foot form	Foot ring	Pk	Cat. No.
Chairs					
Laboratory chair, deformable seat, PU foam, black, tilted 2° forward, with lumbar support adjustable in depth	510 – 645	Castors for hard floors	-	1	132-1104
Laboratory chair, deformable seat, PU foam, black, tilted 2° forward, with high backrest adjustable in depth and lumbar support	510 – 645	Castors with rubber 'tyre'	-	1	132-1300
Laboratory chair, deformable seat, PU foam, black, tilted 2° forward, with lumbar support adjustable in depth	610 – 855	Glides	+	1	132-1105
Laboratory chair, deformable seat, PU foam, black, tilted 2° forward, with high backrest adjustable in depth and lumbar support	610 – 855	Glides	+	1	132-1301
Stools					
Laboratory stool, PU foam, black	510 – 645	Castors for hard floors	-	1	132-1099
Laboratory stool, PU foam, black	610 – 845	Glides	+	1	132-1100

SMART skai chairs



VWR SMART skai chairs with all-inclusive ergonomics. Active, dynamic seating with extraordinary comfort.

- Ergonomically-shaped upholstery guarantees outstandingly comfortable sitting
- Artificial leather upholstery: Washable and disinfectant proof
- Comfortable upholstery with a thickness of 8 cm
- High backrest with integrated lumbar support
- Synchronous technology with weight regulation for active-dynamic sitting
- Height, depth and inclination of seat and backrest individually adjustable
- Polished aluminium 5 star base

GS safety tested.

Description	Height (mm)	Foot form	Foot ring	Pk	Cat. No.
Laboratory chair SMART skai, upholstery black synthetic leather, high backrest, synchronous technology, seat height adjustment	450 – 600	Soft castors	-	1	132-0944
Laboratory chair SMART skai, upholstery blue synthetic leather, high backrest, synchronous technology, seat height adjustment	450 – 600	Soft castors	-	1	132-0945
Laboratory chair SMART skai, upholstery black synthetic leather, high backrest, synchronous technology, seat height adjustment	570 – 830	Glides	+	1	132-0946
Laboratory chair SMART skai, upholstery blue synthetic leather, high backrest, synchronous technology, seat height adjustment	570 – 830	Glides	+	1	132-0947
Laboratory chair SMART skai, upholstery black synthetic leather, high backrest, synchronous technology, seat height adjustment	610 – 870	Stop&Go castors	+	1	132-0948
Laboratory chair SMART skai, upholstery blue synthetic leather, high backrest, synchronous technology, seat height adjustment	610 – 870	Stop&Go castors	+	1	132-0949

Description	Pk	Cat. No.
Accessories		
4-D multifunctional armrest: Height, depth and width adjustable as well as pivoted to fit with SMART and Comfort models	2	132-0702

Chairs, stool, Mr. Lab cleanroom line



132-0421

All chairs and stools are available with a five star base which is made from aluminium.

- Gas lift seat height adjustment
- Adjustable backrest height and inclination
- Sealed upholstery with special foam technology, in artificial leather covering, washable and disinfectant-proof
- Materials suitable for use in cleanrooms
- Smooth, sealed surfaces
- Encapsulated mechanisms
- Electrostatic discharge measure

Glides: The glides are wear-proof plastic and have a large surface area offering a high level of stability.

Castors: Dual-wheel castors with load-sensitive brakes have to be right for the floor. All chairs are supplied with soft castors for hard floors as standard. Hard castors for soft floors are also available as an option.

Stop&Go technology: For safety reasons chairs whose seat is more than 650 mm high must not be movable whilst in use. With innovative Stop&Go castors these chairs are only mobile when no one is sitting on them.

In order to create 'clean' or 'sterile' rooms, measures have to be implemented to prevent products and processes from becoming contaminated. This is a key requirement of the microelectronics sector, the pharmaceuticals industry, microsystem production, optics, medical engineering and healthcare.

VWR cleanrooms chairs meet the requirements of air cleanliness classification 3 in accordance with DIN EN ISO 14644-1, cleanroom classification 1 as defined by US Federal Standard 209E and the specifications of the EU GMP guidelines. Electrostatic discharge protection according to EN 61340-5-1.



132-0424

Description	Height (mm)	Foot form	Foot ring	Pk	Cat. No.
Cleanroom chair, black skai, with permanent contact backrest and aluminium polished base	470 - 620	Castors	-	1	132-0421
Cleanroom chair, black skai, with permanent contact backrest and aluminium polished base	600 - 870	Glides	+	1	132-0422
Cleanroom chair, black skai, with permanent contact backrest and aluminium polished base	640 - 910	Stop&Go castors	+	1	132-0503
Cleanroom stool, black skai, with aluminium polished base	460 - 630	Castors	-	1	132-0424
Description					Cat. No.
Accessories					
Armrest, black					1 PAIR 132-0423

Chairs, Mr. Lab ESD line



132-0419

All chairs and stools are available with a five star base which is made from aluminium.

- Backrest adjustable in height and inclination
- Seat height adjustment by gas-lift
- Aluminium 5 star base
- Upholstery in artificial leather covering, washable and disinfectant-proof

Glides: The glides are wear-proof plastic and have a large surface area offering a high level of stability.

Castors: Dual-wheel castors with load-sensitive brakes have to be right for the floor. All chairs are supplied with soft castors for hard floors as standard. Hard castors for soft floors are also available as an option.

Stop&Go technology: For safety reasons chairs whose seat is more than 650 mm high must not be movable whilst in use. With innovative Stop&Go castors these chairs are only mobile when no one is sitting on them.

ESD (Electro Static Discharge) refers to the electrostatic discharge of charged objects or people. Electrostatic discharge can result in damage to microelectronic components.

Conforms to security level S2, GMP. Soft discharge is achieved on the basis of a fixed and defined discharge resistance of 106 Ω in accordance with DIN EN 61340-5-1.

Description	Height (mm)	Foot form	Foot ring	Pk	Cat. No.
ESD chair made from black skai, with permanent contact backrest and aluminium polished base, plastic parts in black	470 - 620	Castors	-	1	132-0419
ESD chair made from black skai, with permanent contact backrest and aluminium polished base, plastic parts in black	620 - 870	Glides	+	1	132-0420
ESD chair made from black skai, with permanent contact backrest and aluminium polished base, plastic parts in black	640 - 910	Stop&Go castors	+	1	132-0502

Description	Pk	Cat. No.
Accessories		
Armrest, black	1 PAIR	132-0423
Castors	1 SET	132-0427

ESD chairs and stools for cleanrooms



Ergonomic and auto-adaptive chairs or stools made from PU, for cleanrooms

- Auto-adaptive to users and to their positions
- Facilitates an upright posture
- Minimises spinal stress
- Abrasion resistant
- Black steel 5 branch scuff resistant, reinforced base
- Black gas lift column with height adjustment and shock absorber

The innovative seat design consist of a short seat tilted forward, split into 2 parts, based on blue silent-blocks.

Special conductive PU connects to the metal parts of the chairs and stools, allowing electricity to be discharged through carbon castors or pads.

ISO Class 3 Certificate.

Description	Height (mm)	Foot form	Foot ring	Pk	Cat. No.
Chairs					
Cleanroom chair, deformable seat in anti-static black soft PU, tilted 2° forward, with high backrest in PU, adjustable in depth and lumbar support	470 – 605	Heavy-duty anti-static castors	-	1	132-1372
Cleanroom chair, deformable seat in anti-static black soft PU, tilted 2° forward, with lumbar support in PU, adjustable in depth	510 – 645	Anti-static castors	-	1	132-1370
Cleanroom chair, deformable seat in anti-static black soft PU, tilted 2° forward, with high backrest in PU, adjustable in depth and lumbar support	610 – 855	Anti-static glides	-	1	132-1373
Cleanroom chair, deformable seat in anti-static black soft PU, tilted 2° forward, with lumbar support in PU, adjustable in depth	610 – 855	Anti-static glides	+	1	132-1371
Stools					
Cleanroom stool, deformable seat in anti-static black soft PU, tilted 2° forward	510 – 645	Heavy-duty anti-static castors	-	1	132-1368
Cleanroom stool, deformable seat in anti-static black soft PU, tilted 2° forward	610 – 855	Anti-static glides	+	1	132-1369

Chairs, stool, Office chair Smart & Comfort line



132-0686



132-0687



132-0694

A workplace chair with all-inclusive ergonomics: Seat tilt, seat height, backrest height and synchronous mechanism with weight regulation. Synchronous technology with weight regulation: The seat and the backrest tilt are automatically adjusted in line with current posture. There is a continuously variable weight regulation system to ensure that the backrest responds in accordance with how much pressure is applied. It is particularly useful in cases where there are frequent changes in working posture.

- 8 cm thick comfort upholstery
- Integrated lumbar support
- Hard wearing yet soft, comfortable and breathable fabric cover (different colours available)
- Backrest and seat adjustable in height and inclination
- Polished aluminium base

Glides: The glides are wear-proof plastic and have a large surface area for a high level of stability.

Castors: Dual-wheel castors with load-sensitive brakes have to be right for the floor. All chairs are supplied with soft castors for hard floors as standard. Hard castors for soft floors are also available as an option.

Stop&Go technology: For safety reasons chairs whose seat is more than 650 mm high must not be movable whilst in use. With innovative Stop&Go castors these chairs are only mobile when no one is sitting on them.

The Comfort models have a higher backrest than the Smart models.

EN 1335

Description	Height (mm)	Foot form	Foot ring	Pk	Cat. No.
Office chairs, Smart					
Black upholstery	450 - 600	Castors	-	1	132-0686
Blue upholstery	450 - 600	Castors	-	1	132-0687
Dark blue upholstery	450 - 600	Castors	-	1	132-0688
Dark green upholstery	450 - 600	Castors	-	1	132-0689
Red upholstery	450 - 600	Castors	-	1	132-0690
Grey upholstery	450 - 600	Castors	-	1	132-0691
Burgundy upholstery	450 - 600	Castors	-	1	132-0693
Office chairs, Comfort					
Black upholstery	450 - 600	Castors	-	1	132-0694
Blue upholstery	450 - 600	Castors	-	1	132-0695
Dark blue upholstery	450 - 600	Castors	-	1	132-0696
Dark green upholstery	450 - 600	Castors	-	1	132-0697
Red upholstery	450 - 600	Castors	-	1	132-0698
Grey upholstery	450 - 600	Castors	-	1	132-0699
Burgundy upholstery	450 - 600	Castors	-	1	132-0701
Description				Pk	Cat. No.
Accessories					
Armrest, black				1 PAIR	132-0423
4-D multifunctional- armrest: Height, depth and width adjustable as well as pivoted to fit with SMART and Comfort models				2	132-0702
Castors				1 SET	132-0427

Chairs, Lab-Office line with high quality fabric upholstery



These chairs are available with a star base which is made from high strength sectional tubing which is particularly robust and durable. The flat design prevents trip hazards.

- Backrest adjustable in height and inclination
- Seat height adjustment by gas-lift
- Aluminium 5 star base

The perfect combination of comfort and ergonomics for optimal performance. The high quality fabric upholstery is soft, breathable and hard wearing. Anatomically-moulded cushions with pelvis support and integrated lumbar support.

Glides: The glides are wear-proof plastic and have a large surface area offering a high level of stability.

Castors: Dual-wheel castors with load-sensitive brakes have to be right for the floor. All chairs are supplied with soft castors for hard floors as standard. Hard castors for soft floors are also available as an option.

Stop&Go technology: For safety reasons chairs whose seat is more than 650 mm high must not be movable whilst in use. With innovative Stop&Go castors these chairs are only mobile when no one is sitting on them.



Description	Height (mm)	Foot form	Foot ring	Pk	Cat. No.
Office chair, made from blue fabric	470 - 620	Castors	-	1	132-0418
Description				Pk	Cat. No.
Accessories					
Armrest, black				1 PAIR	132-0423
Castors				1 SET	132-0427

Footrest, relax



- Adjustable in height and inclination
- Recess for foot switch
- The ribbed platform ensures a firm grip and avoids slipping

In compliance with DIN standard 4556, GS tested.

Description	Colour	WxDxH (mm)	Weight (kg)	Pk	Cat. No.
Footrest, relax	Light grey	490x350x190	2,4	1	139-9923



Footrest, Trittboy

Made from steel tubing, with nickel-plated finish.

- Four adjustable height levels
- Stable and non slip

Colour: Black tread

Minimum height: 60 mm - maximum height: 145 mm



Description	Colour	WxDxH (mm)	Pk	Cat. No.
Footrest, Trittboy	Black	345x335x16	1	139-9922

Footrest, swing



A stylish footrest that is fully adjustable from a normal seated position. Ideal for use at shared workstations.

- Ribbed platform provides extra grip

Description	Colour	WxDxH (mm)	Weight (kg)	Pk	Cat. No.
Footrest, swing	Black	4770x320x105	1,7	1	139-9805



PCR Workstation



The VWR® PCR Workstation offers twin decontamination action by UV inactivation of airborne and surface-bound contaminants and, therefore, represents an ideal environment for PCR sample preparation and other sensitive protocols. Housing in stainless steel, front and side panels made of 8 mm Makrolon® protects from UV irradiation.

- Active decontamination of work surface during non working time by UV irradiation
- Additional inactivation of aerosol-bound contaminants by shielded UV Air Recirculator during operation
- Contaminant prevention thanks to antimicrobial stainless steel work surface
- Function indicator for UV Air Recirculator tube
- Displays operating time of UV tubes to enable well-timed replacement for constant UV intensity
- Removable shelves provide additional storage space for reaction tubes, pipettes or racks at the rear panel
- Four internal power outlets for operating lab equipment such as mini centrifuges or vortexers allow the combination of several working steps without interruption of the workflow
- Electromagnetic safety mechanism stops UV irradiation if front panel is opened

Large working area for stress-free working: Providing ample space and the possibility to accommodate bench top equipment, the VWR® PCR Workstation allows the combination of several working steps without change of location, thereby minimising the risk of cross-contamination.

TÜV-tested, safety certified. Made in Germany.



Model	PCR Workstation		
Lighting	2 UV tubes internal (254 nm, 25 W each) 1 UV tube in UV Air Recirculator (254 nm, 8 W) 1 white light tube internal (15 W)		
Shelves	Removable		

Description	WxDxH ext. (mm)	Pk	Cat. No.
VWR® PCR Workstation, work surface (WxD): 720x540 mm, EU-plug	750x620x780	1	732-2541
VWR® PCR Workstation, work surface (WxD): 720x540 mm, UK-plug	750x620x780	1	732-2542

Description	Pk	Cat. No.
Accessories		
Light source, interior, white light, 15 W	1	732-2543
Light source, interior, UV light (254 nm), 15 W	2	732-2544
Light source, air recirculator, UV light (254 nm), 8 W	1	732-2545
Dust filter for air recirculator (8 W)	10	732-2546

Weighing table, granite, for micro and analytical balances

High performance microbalances and analytical balances should be placed in an environment where interfering vibration is at a minimum, to give fast stabilisation of the reading. The VWR weighing table features a sturdy, stable frame.

- Granite weighing surface
- Adjustable feet
- Meets highest anti-vibration standards, with four vibration dampers
- Facilitates fast stabilisation times

Installation not included.

Description	Colour	WxDxH (mm)	Weight (kg)	Pk	Cat. No.
Balance table, granite, with vibration dampers and adjustable feet	Grey	900x600x760	267	1	611-2549



Anti-vibration weighing surface

Polished granite

Four absorbing feet dampen the transmission of extraneous vibrations to this slab making it suitable, not only for high precision weighing, but also microscopy and other similar delicate instrumentation.

Description	WxDxH (mm)	Weight (kg)	Pk	Cat. No.
Anti-vibration weighing surface, granite	560x450x65	44	1	611-2536



LED illuminated bench magnifiers



These illuminated magnifiers can be fixed on a table or on a bench. The square magnifier has 80 LED lights, while the round models has 90.

- Modern design
- Articulated arm (800 mm) for a 200° rotation and better manipulation
- Cover to protect the lens from dust and scratches when the magnifier is not being used

Includes a standard table bracket (WxDxH: 30x40x64 mm).

Description	Magnification	Pk	Cat. No.
Square magnifier, EU-plug	1,5x	1	630-1947
Square magnifier, UK-plug	1,5x	1	630-1948
Round magnifier, EU-plug	1,5x	1	630-1950
Round magnifier, UK-plug	1,5x	1	630-1951

Description	Pk	Cat. No.
Accessories		
Replacement LED for 80 LED magnifiers	1	630-1962
Replacement LED for 90 LED magnifiers	1	630-1963

LED illuminated table magnifiers



These illuminated desk top magnifiers - gooseneck or articulated arm type - include highly polished glass.

- Large viewing area
- Cover protects lens from dust and scratches
- Articulated arm models provide 200° head rotation

Description	Magnification	Pk	Cat. No.
Gooseneck magnifier, 60 LED, EU-plug	1,75x ; 3x	1	630-1953
Gooseneck magnifier, 60 LED, UK-plug	1,75x ; 3x	1	630-1954
Rectangular base magnifier, 60 LED, EU-plug	1,75x ; 3x	1	630-1956
Rectangular base magnifier, 60 LED, UK-plug	1,75x ; 3x	1	630-1957
Round base magnifier, 40 LED, EU-plug	1,75x ; 3x	1	630-1959
Round base magnifier, 40 LED, UK-plug	1,75x ; 3x	1	630-1960

Description	Pk	Cat. No.
Accessories		
40 LED illumination	1	630-1965
60 LED illumination	1	630-1964

Stepladder, mobile



Stepladder, mobile in stainless steel, with moon castors. While weight is applied on the first step, wheels retract and the moon cups with cleanroom approved nylon guards sit firmly on the floor, ensuring a safe climb.

- All the ladders can be manoeuvred by a single person
- Ideal if they have to be moved often

Height of top step from ground level when weight is applied are as follows:

1 step ladder-201 mm

2 step ladder-520 mm

3 step ladder-609 mm

4 step ladder-782 mm

Description	Height (mm)	WxDxH (mm)	Pk	Cat. No.
Mobile step ladder, 1 step, with retractable castors	200	600x300x200	1	139-0029
Mobile step ladder, 2 step, with retractable castors	600	805x500x600	1	139-0030
Mobile step ladder, 3 step, with retractable castors	1509	632x832x1509	1	139-0031

Folding steps



Foldable steps with handle, sturdy double side-rails, grooved aluminium steps (205x360 mm).

- Easy fold away
- Non slip and robust for a safe working
- Measures only 85 mm in depth when folded

GS tested; in accordance with European standard DIN EN 14183.

Action manual BGI 694 and the applicable accident prevention regulations.

Description	No. of steps	Working reach up to approx.	Base spread	Height of the top step	Outside width	Weight (kg)	Pk	Cat. No.
Folding step with 3 rungs and 1 handle	3	2,70 m	0,82 m	0,64 m	570 mm	7	1	139-0770



Folding steps



Foldable steps with handle, sturdy double side rails, grooved aluminium steps (205x360 mm), foldable bracket and tray.

- Non slip feet
- Robust and safe
- Only 125 mm in depth when folded
- 15 year quality guarantee

GS tested.



Description	No. of steps	Working reach up to approx.	Base spread	Height of the top step	Weight (kg)	Pk	Cat. No.
Folding step	3	2700 mm	820 mm	640 mm	9	1	139-0426



Single-sided step ladders, with castors

Ascendable on one side, grooved steps with a depth of 80 mm with durable step-to-side-rail connection, non slip nivello® ladder feet. High standing safety ensured by conical construction.

- Connection on both sides by two detachable metal struts
- Four self-locking fixed castors at the side rail ends
- Spacing between steps: 235 mm
- Ladder inclination: 20°
- Non slip nivello® ladder shoes
- Ergopad grip zone
- Safe cap joint
- Durable step to side rail connection

GS tested; in accordance with European standard DIN EN 14183.

Action manual BGI 694 and the applicable accident prevention regulations.

Description	No. of steps	Working reach up to approx.	Height (mm)	Platform height	Weight (kg)	Pk	Cat. No.
Step ladder with 3 steps	3	2750 mm	1370	700 mm	5,3	1	139-0769
Step ladder with 3 steps	3	2750 mm	1370	700 mm	8,1	1	139-0650
Step ladder with 3 steps and step-pads	3	2750 mm	1370	700 mm	8,1	1	139-0652
Step ladder with 3 steps and slip resistant clip-step inserts	3	2750 mm	1370	700 mm	8,1	1	139-0651

Laboratory step ladders, aluminium

Step ladders made from aluminium. 15 years of quality guaranteed.

- Slim-line construction for space saving storage
- High stability due to conical design
- Payload up to 150 kg
- Corrugated aluminium steps (80 mm deep), ensuring a safe stance
- Large, corrugated platform for safe working
- Non slip corrugated ladder feet to ensure full contact with the underlying surface
- Extra light and easy to use

Colour	No. of steps	Working reach up to approx.	Weight (kg)	Pk	Cat. No.
Aluminium	3	2,75 m	4,7	1	139-0469
Aluminium	4	3 m	5,8	1	139-0470
Aluminium	5	3,25 m	6,9	1	139-0471
Aluminium	6	3,50 m	8	1	139-0472
Aluminium	7	3,75 m	9,4	1	139-0473



Transport trolleys

Transport trolley with high quality stainless steel construction, complete with stable tubes of 25x25 mm, and two smooth-welded shelves.

- Zinc-chrome plate steel housing, anti-static with bumpers
- Chassis complete with four swivel castors, front two lockable

Description	WxDxH (mm)	Shelves	Pk	Cat. No.
Multi-purpose trolley, with 4 swivel castors (2 with brake)	625x425x800	2	1	139-0771
Multi-purpose trolley, with 4 swivel castors (2 with brake)	1200x625x800	2	1	139-0772



Transport trolleys

Chrome plated steel frame with HDPE trough shaped trays.

- Stable, tough and versatile, easy to clean
- Load capacity: Max. 65 kg for each tray; combined 180 kg
- Rubber wheels to protect floors

The troughs are 70 mm deep, enabling liquids to be transported safely retaining any spillage. The height of the middle tray can be adjusted. Available in two sizes, with two or three trays.

WxDxH (mm)	Shelves	Colour	Pk	Cat. No.
Shelf size: 830x545 mm				
830x545x1040	2	Blue	1	139-0444
830x545x1040	2	Grey	1	139-9991
Shelf size: 985x690 mm				
985x690x1040	2	Grey	1	139-9993
985x690x1040	3	Grey	1	139-9994

Utility carts, MW series

These steel or chrome-plated carts consist of Super Erecta shelves with plastic split sleeves, two one-piece handles of matching finish, and designated castors with donut bumpers. These multipurpose trolleys are perfect for carrying out an enormous variety of tasks.

- Ergonomically-designed for low effort transport
- Precise handling to manoeuvre in close quarters
- Sturdy aluminum corner castings firmly secure solid shelves to posts
- Stable, tough, simple and easy to clean
- MW carts use light-duty resilient castors



Description	WxDxH (mm)	Weight (kg)	Shelves	Pk	Cat. No.
Utility carts with solid shelves (stainless steel)					
MW103	610x457x965	18,1	2	1	139-0776
MW104	760x457x965	20,4	2	1	139-0777
MW105	914x457x965	21,8	2	1	139-0778
MW106	914x530x990	24,5	2	1	139-0779
MW108	914x610x990	27,2	2	1	139-0780
MW203	610x457x965	18,1	3	1	139-0372
MW204	760x457x965	20,4	3	1	139-0373
MW205	914x457x965	21,8	3	1	139-0374
MW206	914x530x990	24,5	3	1	139-0375
MW208	914x610x990	27,2	3	1	139-0781
Utility carts with 1 solid shelf (stainless steel) and 2 wire shelves (chrome plated)					
MW401	610x457x965	19,5	1+2	1	139-0077
MW402	760x457x965	21,3	1+2	1	139-0078
MW403	914x457x965	23,6	1+2	1	139-0079
MW404	914x530x990	26,8	1+2	1	139-0080
MW406	914x610x990	29,5	1+2	1	139-0081
Utility cart with 2 wire shelves (stainless steel)					
MW602	610x457x965	15,4	2	1	139-0376
MW604	760x457x965	16,8	2	1	139-0377
MW606	914x457x965	17,7	2	1	139-0378
MW608	914x530x990	20	2	1	139-0379
MW612	914x610x990	20,9	2	1	139-0380
Utility cart with 2 wire shelves (chrome plated)					
MW601	610x457x965	15,4	2	1	139-0381
MW603	760x457x965	16,8	2	1	139-0382
MW605	914x457x965	18,1	2	1	139-0383
MW607	914x530x990	20	2	1	139-0384
MW611	914x610x990	21,3	2	1	139-0385

Description	WxDxH (mm)	Weight (kg)	Shelves	Pk	Cat. No.
Utility carts with 3 wire shelves (stainless steel)					
MW702	610x457x965	17,2	3	1	139-0388
MW704	760x457x965	19,5	3	1	139-0389
MW706	914x457x965	20,4	3	1	139-0390
MW708	914x530x990	22,7	3	1	139-0391
MW712	914x610x990	25,9	3	1	139-0392
Utility carts with 3 wire shelves (chrome plated)					
MW701	610x457x965	17,6	3	1	139-0386
MW703	760x457x965	19,5	3	1	139-9995
MW705	914x457x965	21,3	3	1	139-9996
MW707	914x530x990	22,7	3	1	139-0387
MW711	914x610x990	27,7	3	1	139-9997



Transport trolleys, myCart™ series

myCart™ series utility carts are made from corrosion proof polymer shelves and are suitable for daily use.

- Smooth surfaces easily wipe clean and do not stain
- Additional leg room provides more distance between the shelves and the user
- Available in two- and three- shelf configurations
- 11 mm deep ship's edge lip around shelf perimeter

Description	WxDxH (mm)	Weight (kg)	Shelves	Colour	Pk	Cat. No.
Transport trolley, base model, chrome posts, with handle	799x465x899	12,5	2	Black	1 SET	139-0782
Transport trolley, base model, chrome posts, with handle, with antimicrobial protection	799x465x899	12,5	2	Blue	1 SET	139-0783
Transport trolley, base model, chrome posts, with handle	799x465x899	12,5	2	Grey	1 SET	139-0785
Transport trolley, base model, chrome posts, with handle	799x465x899	15	3	Black	1 SET	139-0786
Transport trolley, base model, chrome posts, with handle, with antimicrobial protection	799x465x899	15	3	Blue	1 SET	139-0787
Transport trolley, base model, chrome posts, with handle	799x465x899	15	3	Grey	1 SET	139-0788
Transport trolley, base model, chrome posts, with handle	873x597x899	15,7	2	Black	1 SET	139-0789
Transport trolley, base model, chrome posts, with handle, with antimicrobial protection	873x597x899	15,7	2	Blue	1 SET	139-0430
Transport trolley, base model, chrome posts, with handle	873x597x899	15,7	2	Grey	1 SET	139-0429
Transport trolley, base model, chrome posts, with handle	873x597x899	19,5	3	Black	1 SET	139-0790
Transport trolley, base model, chrome posts, with handle, with antimicrobial protection	873x597x899	19,5	3	Blue	1 SET	139-0432
Transport trolley, base model, chrome posts, with handle	873x597x899	19,5	3	Grey	1 SET	139-0431
Transport trolley, base model, chrome posts, with handle	1022x703x930	18,4	2	Black	1 SET	139-0791
Transport trolley, base model, chrome posts, with handle, with antimicrobial protection	1022x703x930	18,4	2	Blue	1 SET	139-0433
Transport trolley, base model, chrome posts, with handle	1022x703x930	18,4	2	Grey	1 SET	139-0434
Transport trolley, base model, chrome posts, with handle	1022x703x930	23,6	3	Black	1 SET	139-0792
Transport trolley, base model, chrome posts, with handle, with antimicrobial protection	1022x703x930	23,6	3	Blue	1 SET	139-0435
Transport trolley, base model, chrome posts, with handle	1022x703x930	23,6	3	Grey	1 SET	139-0436

Utility carts, SP series

These carts are designed to transport heavy loads. They consist of Super Erecta wire shelves with plastic split sleeves, two one-piece handles with matching finish and designated castors with donut bumpers. Made from steel, Super Erecta brite finish.

- Two-tier model
- 5M resilient rubber castors, castors have plated horns and axles

Description	WxDxH (mm)	Weight (kg)	Shelves	Pk	Cat. No.
Utility cart, Super Erecta brite finish, with 5M castors	1219x610x990	24	2 wire shelves	1	139-0773



Utility carts



These cleanroom carts are designed to withstand heavy loads.

- Up to 60 kg per shelf
- Integral raised ledges prevent parts from rolling off and are easy to wipe down
- Choice of 2 or 3 shelves
- Solid or perforated
- White nylon 130 mm (5») cleanroom castors (2 with brakes)

Stainless steel 304 with electropolished finish

Shipped flat-packed with assembly instructions.

WxDxH (mm)	Shelves	Pk	Cat. No.
Cleanroom carts with solid shelves			
550x550x990	2	1	139-0588
550x550x990	3	1	139-0793
865x550x990	2	1	139-0589
865x550x990	3	1	139-0590
Cleanroom carts with perforated shelves			
550x550x990	3	1	139-0794
865x550x990	2	1	139-0591
865x550x990	3	1	139-0795

Transport trolleys for glassware, MetroMax i®

The MetroMax i® glassware carts are available in two and four basket capacity and are ideal for transporting glassware safely in the laboratory.

- The solid bottom shelf provides spill retention and can be removed from shelf frame for easy cleaning

Castors must be ordered separately.

Description	WxDxH (mm)	Weight (kg)	Pk	Cat. No.
Cart, 2 basket unit	914x457x1016	17	1	139-0062
Cart, 4 basket unit	914x457x1524	37	1	139-0063
Cart, 2 tote kit unit	914x457x1016	16,3	1	139-0069
Cart, 4 tote kit unit	914x457x1524	46	1	139-0070



Description	Pk	Cat. No.
Polymer stem castors		
Polymer stem castors, brake, polyurethane	1	139-0802
Polymer stem castors, swivel, polyurethane (antimicrobial product)	1	139-0803
Polymer stem castors, brake, polyurethane (antimicrobial product)	1	139-0050
Standard stem castors		
Standard stem castors, stem/swivel, polyurethane	1	139-0047
Standard stem castors, stem/brake, polyurethane	1	139-0046
Standard stem castors, stem/swivel, resilient donut	1	139-0801
Standard stem castors, stem/brake, resilient donut	1	139-0045
Stainless steel cart-washable stem castors		
Stainless steel cart-washable stem castors, stem/swivel, polyurethane	1	139-0048
Stainless steel cart-washable stem castors, stem/brake, polyurethane	1	139-0083
Accessories		
Cart cover, for short unit, clear, 914x457x889 mm	1	139-0064
Cart cover, for tall unit, clear, 914x457x1320 mm	1	139-0065
Cart cover, for short unit, blue, 914x457x889 mm	1	139-0067
Cart cover, for tall unit, blue, 914x457x1320 mm	1	139-0068
Vinyl coated baskets, 2, 901x457x178 mm	1	139-0061
Tote glassware transport kit	1	139-0066

Transport trolleys



Transport trolley with high quality stainless steel construction, complete with 25 mm deep shelves. Shelves have bevelled edges for user comfort.

- Bumper-type wall deflectors
- Anti-static tyres
- Chassis complete with four swivel castors, front two lockable



Description	WxDxH (mm)	Shelves	Pk	Cat. No.
Transport trolley MSW 6x4/2, standard version, 2 shelves WxD: 600x400 mm	700x500x895	2	1	139-1025
Transport trolley MSW 8x5/2, standard version, 2 shelves WxD: 800x500 mm	900x600x960	2	1	139-9940
Transport trolley MSW 10x6/2, standard version, 2 shelves WxD: 1000x600 mm	1095x695x960	2	1	139-1022
Transport trolley MSW 8x5/2, heavy duty, 2 shelves WxD: 800x500 mm	900x600x1030	2	1	139-1024
Transport trolley MSW 8x5/3, standard version, 3 shelves WxD: 800x500 mm	900x600x960	3	1	139-9909
Transport trolley MSW 10x6/3, standard version, 3 shelves WxD: 1000x600 mm	1095x695x960	3	1	139-1023
Transport trolley MSW 8x5/4, standard version, 4 shelves WxD: 800x500 mm	900x600x1280	4	1	139-9945



Transport trolleys

These trolleys made of stainless steel are mobile on four swivel castors in galvanised chromed housings, two of them with total brake, Ø 125 mm.

Platform trolley: Flat metal platform, folded down by 80 mm all around with two plastic bumpers at the front. It has four swivel castors, two of them with total brake, Ø 125 mm

Single-ended trolley: Platform manufactured from round tube, Ø 25 mm, with welded-on stainless steel plate. It has two swivel castors with total brake, two fixed castors, Ø 200 mm

Description	WxDxH (mm)	Weight (kg)	Pk	Cat. No.
Platform trolley				
Platform trolley PW, load capacity 120 kg	723x605x1022	12,2	1	139-9904
Platform trolley PW	983x605x1020	17,1	1	139-9906
Single ended trolley				
Single ended trolley STWP	970x500x940	42	1	139-9902
Single ended trolley STWP	1120x700x940	52	1	139-9903



Workstation dispensers



The acrylic combination accessory dispenser holds the most commonly used laboratory disposables in one central location. Dispenser stores tubes, tips, gloves, Kimwipes®, and Parafilm®.

- Two compartments with hinged lid
- One compartment for wipe box and one compartment for tape
- One large compartment for glove or wipe box

Description	WxDxH (mm)	Colour	Pk	Cat. No.
Workstation dispenser	265 x 137 x 440	Clear	1	216-2006

Disposal box for pipette tips



This acrylic box is ideal for the one-handed disposal of used pipette tips. The 10 mm thick acrylic provides protection from β radiation.

- Quick flip of port cover makes tip ejection easy
- Eliminates any need to handle pipette tips used with hazardous fluids
- Large hinged lid allows easy disposal of contents

WxDxH: 154x171x206 mm

Description	Pk	Cat. No.
Disposal box for pipette tips, clear acrylic, with flip cover, 154x171x206 mm	1	212-1833

Dispensing bins



Acrylic

Convenient, space saving bins are ergonomically designed with large angled openings to facilitate the removal of items. Bins feature non skid rubber feet for added bench top stability as well as pre-drilled holes to allow for wall mounting (hardware included).

- Lid tilts back into locking position to remain open for repeated use
- Lid features an oversized lip to provide easy access, even to gloved hands
- Green bin is supplied with a removable insert divider to allow customisation of the interior (one or two compartments)

Description	Colour	WxDxH (mm)	Pk	Cat. No.
Dispensing bin, 1 compartment	Blue	115x155x170	1	631-0720
Dispensing bin, 2 compartments	Green	180x155x165	1	631-0721

Drawer organiser



PS.

- One wide compartment: 430 mm long, 110 mm wide
- Two compartments each 430 mm long, 80 mm wide
- Smooth surfaces with rounded corners for easy cleaning

For well organised storage of various different items in drawers, on shelves or at the workstation. Three long compartments for pipettes, thermometers and spatulas.

WxDxH (mm)	Colour	Pk	Cat. No.
460x307x63	White	1	631-0727

Workstation dispensers



Crystal-clear bins for small items inventory, like pipette tips, tubes, fittings, finger cots, and other bulk small items.

- Flip-top covers on generous openings make filling and dispensing easy and keep products dust-free
- 4 sizes available to cater to a wide variety of storage needs
- Includes non-skid rubber feet

Description	WxDxH (mm)	Pk	Cat. No.
Workstation dispensing bin, one compartment	180x175x235	1	216-2002
Workstation dispensing bin, two compartments	139x160x368	1	216-2003
Workstation dispensing bin, three compartments	210x155x212	1	216-2004
Workstation dispensing bin, four compartments	278x158x217	1	216-2005

Storage boxes




PS, clear.

- Stackable, transparent boxes
- For storage and presentation
- Six different sizes

These boxes are an ideal choice for laboratories due to their very high transparency and excellent durability.

WxDxH (mm)	Pk	Cat. No.
60x45x50	1	216-4332
90x60x50	1	216-4333
120x90x50	1	216-4334
180x120x75	1	216-4000
240x180x100	1	216-4002
360x240x140	1	216-4003

Tilt storage boxes



ABS


Each element contains a number of removable transparent drawer suitable for a wide range of items.

- Non toxic, suitable for use with foodstuffs
- Strong outer body high impact polystyrene plastic and inner compartments
- Can be fixed directly to the wall



Description	WxDxH (mm)	Colour	Pk	Cat. No.
1 drawer	122x133x160	White	4	216-4361
2 drawers	600x299x353	Beige	1	216-4355
3 drawers	600x197x240	Beige	1	216-4348
4 drawers	600x168x207	Beige	1	216-4356
5 drawers	600x133x164	Beige	1	216-4357
6 drawers	600x91x113	Beige	1	216-4358
9 drawers	600x62x77	Beige	1	216-4359

Drawers, transparent



Longitudinal, transverse and semi transverse

System of modular drawer cabinets made from PS. For well organised storage of many items, on shelves or at a workstation.

- Total visibility in the drawer
- Each drawer can be assembled horizontally and vertically
- Each drawer can be divided with separators

Capacity (l)	WxDxH (mm)	Colour	Pk	Cat. No.
1,9	235x168x72	Transparent	1	216-4346
2	275x166x61	Transparent	1	216-4347
7	365x210x125	Transparent	1	216-4349
10	410x265x125	Transparent	1	216-4326
30	542x322x216	Transparent	1	216-4370

Description	Pk	Cat. No.
Partitions for storage boxes		
Longitudinal separator for drawer 216-4346 (1901)	1	216-4327
Transverse separator for drawer 216-4346 (1902)	1	216-4328
Semi-transverse separator for drawer 216-4346 (1903)	1	216-4329
Longitudinal separator for drawer 216-4347 (2001)	1	216-4336
Transverse separator for drawer 216-4347 (2002)	1	216-4337
Semi-transverse separator for drawer 216-4347 (2003)	1	216-4338
Longitudinal separator for drawer 216-4349 (7001)	1	216-4339
Transverse separator for drawer 216-4349 (7002)	1	216-4340
Semi-transverse separator for drawer 216-4349 (7003)	1	216-4341
Longitudinal separator for drawer 216-4326 (10001)	1	216-4342
Transverse separator for drawer 216-4326 (10002)	1	216-4343
Semi-transverse separator for drawer 216-4326 (10003)	1	216-4344

Glove box holder



Glove boxes holder in durable resin to resist to harsh cleaning solutions.

- Wall-mountable design for easy use
- Modular design for horizontal stacking of multiple boxes
- Universal fit for many of the major brands

Suitable for use in cleanrooms, hospitals and microbiological working areas.

Mounting screws included.

Description	WxDxH (mm)	Colour	Pk	Cat. No.
Glove boxes holder	254x156x97	White	3	611-3119

Glove box holders



Holders made from tough clear acrylic. You can clearly see the size and material of the gloves within.

- Can hold one, two or three glove boxes
- Polyester foam insert for a secure storage
- Strong, attractive, scratch resistant material

Wall-mount screws included.

Description	WxDxH (mm)	Colour	Pk	Cat. No.
Glove box (1) holder	141x110x260	Clear	1	113-8041
Glove boxes (2) holder	295x110x260	Clear	1	113-8042
Glove boxes (3) holder	295x110x385	Clear	1	113-8043

Glove box holders



Epoxy coated steel wire racks.

- Available as a single box or three-box holder
- Can hold most major brands of boxes
- Easy mounting on the wall

Screws are included.

Description	WxDxH (mm)	Pk	Cat. No.
Glove box holder	140x108x208	1	211-0343
Glove boxes (3) holder	250x108x454	1	211-0344

Glove box holders

These wall-mountable glove box holders are made from polystyrene with Microban® additive.

- Microban® additive to prevent growth of moulds and bacteria
- High impact polystyrene, resistant to harsh cleaning solutions
- Modular design allows stacking to save space

Suitable for use in cleanrooms, hospitals and microbiological work areas.

Mounting screws included.



Description	WxDxH (mm)	Colour	Pk	Cat. No.
Glove box holder, Microban®	254x156x97	Blue	3	211-0371

Wall mount for four glove or clothing boxes

The adjustable holding points provide true flexibility to hold up to four different boxed protective products and other lab essentials. Ideal for different boxed sizes of gloves, face masks, ear plugs or shoe covers.

- Four holding positions can hold up to four boxes in one holder
- Wall or door-mountable, mounting parts included
- No restrictions regarding box size: Open ended holding points are able to hold longer and wider boxes
- Easy to clean, with low profile



Type	LxW (mm)	Colour	Pk	Cat. No.
Wall holder, adjustable, up to 4 positions	762x152	White	1	113-8040

designed for
every day use



VWR® LAB BRIC-A-BRAC


Never know where to put anything?
Need extra storage?
New batteries?

All the things you didn't know you needed until you couldn't find them or didn't have them.

Make lab life easier with VWR.

For your copy visit vwr.com

Mini equipment log books



'FLEX' covers with coating, for durability in laboratory environments. Acid-free, archive-safe paper (90 gsm).

- Pages specifically designed for logging and recording equipment data
- Tough coated cover

Format: A5 (14,8x21 cm) or A6 (10,5x14,8 cm), 96 pages, numbered.

All books made in accordance with GLP, GMP and GCP guidelines.

Description	Colour	Pk	Cat. No.
Mini equipment log book, A5, 96 pages	Red	1	818-0108
Mini equipment log book, A6, 96 pages	Red	1	818-0109

Log book



Professional grade log book for tamper-evident record keeping.

- Blue cover
- 48 pages

Description	Colour	Pk	Cat. No.
Log book, A4, 48 pages	Blue	1	818-9994

Mini laboratory notebooks


Imitation leather for use in laboratory environments. Acid-free archive safe paper (90 gsm). Includes laboratory notebook guidelines page, table of contents and user data page. English language only.

- Scientific 5 mm grid or ruled format with date and witness signature blocks
- 'FLEX' cover with coating, for durability

Format: A5 (14,8x21 cm) or A6 (10,5x14,8 cm), 96 pages

Colour: Blue

All books made in accordance with GLP, GMP and GCP guidelines.




Description	Colour	Type	Pk	Cat. No.
Mini lab book, A5	Blue	Grid	1	818-0100
Mini lab book, A6	Blue	Grid	1	818-0101
Mini lab book, A5	Blue	Ruled	1	818-0110
Mini lab book, A6	Blue	Ruled	1	818-0111

Laboratory log books, special purpose


Imitation leather cover with coating, for durability in laboratory environments. Acid-free, archive-safe paper (90 gsm). Used to record and log data relating to laboratory equipment, temperature, balance calibration, reagents, animal maintenance and sample receipts.

- All pages for logging documentation have specific column headings for easy data entry
- Pages sewn in sections, books have round binding so that they lay flat
- Waterproof and tough

All books made in accordance with GLP, GMP and GCP guidelines.



Description	Colour	Pk	Cat. No.
Log book, laboratory equipment, 96 pages	Red	1	818-0102
Log book, animal maintenance, 168 pages	Black	1	818-0103
Log book, balance calibration, 168 pages	Green	1	818-0104
Log book, reagents, 168 pages	Dark red	1	818-0105
Log book, sample receipt, 168 pages	Red	1	818-0106
Log book, temperature, 168 pages	White	1	818-0107



Laboratory notebooks

Imitation leather cover with coating, for durability in laboratory environments. Acid-free archive safe paper (90 gsm). Includes laboratory notebook guidelines page, table of contents and user data page. Available in different language formats: English only, English-French, English-German, English-Spanish and English-Italian.

- Pages are sewn in sections with round binding so that book lays flat
- Scientific 5 mm grid or ruled format with date and witness signature blocks
- Waterproof

Format: A4 (21×29,7 cm) or A5 (14,8×21 cm). A4 books available in three sizes: 72, 120 and 216 pages. A5 available with 96 pages.

Colour: Blue

All books made in accordance with GLP, GMP and GCP guidelines.

Description	Colour	Type	Version	Pk	Cat. No.
Lab book, A4, 72 pages	Blue	Grid	English only	1	818-0068
Lab book, A4, 120 pages	Blue	Grid	English only	1	818-0069
Lab book, A4, 216 pages	Blue	Grid	English only	1	818-0070
Lab book, A5, 96 pages	Blue	Grid	English only	1	818-0071
Lab book, A4, 72 pages	Blue	Ruled	English only	1	818-0072
Lab book, A4, 120 pages	Blue	Ruled	English only	1	818-0073
Lab book, A4, 216 pages	Blue	Ruled	English only	1	818-0074
Lab book, A5, 96 pages	Blue	Ruled	English only	1	818-0075
Lab book, A4, 72 pages	Blue	Grid	English-French	1	818-0076
Lab book, A4, 120 pages	Blue	Grid	English-French	1	818-0077
Lab book, A4, 216 pages	Blue	Grid	English-French	1	818-0078
Lab book, A5, 96 pages	Blue	Grid	English-French	1	818-0079
Lab book, A4, 72 pages	Blue	Ruled	English-French	1	818-0080
Lab book, A4, 120 pages	Blue	Ruled	English-French	1	818-0081
Lab book, A4, 216 pages	Blue	Ruled	English-French	1	818-0082
Lab book, A5, 96 pages	Blue	Ruled	English-French	1	818-0083
Lab book, A4, 72 pages	Blue	Grid	English-German	1	818-0084
Lab book, A4, 120 pages	Blue	Grid	English-German	1	818-0085
Lab book, A4, 216 pages	Blue	Grid	English-German	1	818-0086

Description	Colour	Type	Version	Pk	Cat. No.
Lab book, A5, 96 pages	Blue	Grid	English-German	1	818-0087
Lab book, A4, 72 pages	Blue	Ruled	English-German	1	818-0088
Lab book, A4, 120 pages	Blue	Ruled	English-German	1	818-0089
Lab book, A4, 216 pages	Blue	Ruled	English-German	1	818-0090
Lab book, A5, 96 pages	Blue	Ruled	English-German	1	818-0091
Lab book, A4, 72 pages	Blue	Grid	English-Spanish	1	818-0092
Lab book, A4, 120 pages	Blue	Grid	English-Spanish	1	818-0093
Lab book, A4, 216 pages	Blue	Grid	English-Spanish	1	818-0094
Lab book, A5, 96 pages	Blue	Grid	English-Spanish	1	818-0095
Lab book, A4, 72 pages	Blue	Ruled	English-Spanish	1	818-0096
Lab book, A4, 120 pages	Blue	Ruled	English-Spanish	1	818-0097
Lab book, A4, 216 pages	Blue	Ruled	English-Spanish	1	818-0098
Lab book, A5, 96 pages	Blue	Ruled	English-Spanish	1	818-0099
Lab book (21,6×14,0 cm), 120 pages	Blue	Grid	English-Italian	1	818-0112
Lab book, A4, 216 pages	Blue	Grid	English-Italian	1	818-0113
Lab book, A4, 120 pages	Blue	Ruled	English-Italian	1	818-0114
Lab book (21,6×14,0 cm), 216 pages	Blue	Ruled	English-Italian	1	818-0115
Lab book, XL (22,5×34,5 cm), 144 pages	Blue	Grid	EN/FR/IT/SP/DE	1	818-0116
Lab book, XL (22,5×34,5 cm), 144 pages	Blue	Ruled	EN/FR/IT/SP/DE	1	818-0117

Laboratory markers



Versatile lab markers are designed to write on wet or dry surfaces of disposable plastic material, glassware or porcelain.

- Ink will not erase when it comes in contact with water or alcohol
- Ink dries instantly, allowing immediate handling
- Can be removed with solvents (except alcohol)

Description	Colour	Pk	Cat. No.
Lab marker, fine tip	Black	10	811-0065
Lab marker, fine tip	Blue	10	811-0066
Lab marker, fine tip	Red	10	811-0067
Lab marker, fine tip	Green	10	811-0068

Permanent marker pens



Perfect for slides, sample bags and glassware.

- Write on glass, plastic, porcelain, film, rubber, wood and paper
- Permanent ink dries in seconds

Assorted colours: 6 black, 2 red, 2 blue, 2 green

Description	Colour	Pk	Cat. No.
Marker pen	Assorted	12	811-0055
Marker pen	Black	6	811-0056
Marker pen	Red	6	811-0057
Marker pen	Blue	6	811-0058
Marker pen	Green	6	811-0059



Wipes, pre-saturated

These non woven wipes are constructed of cellulose and polyester and moistened with a blend of pure, reagent-grade 70% isopropyl alcohol and 30% reagent-grade deionised water. They are ideal for wiping optical parts, delicate glassware, electrodes, cuvettes, microscopes, or lenses.


- Suitable for all applications in environmentally-sensitive areas
- Very strong, soft and lint-free
- Will not scratch surfaces
- Individually dispensed wipes are handy and ready for instant use
- Reclosable bench top container ensures a clean wipe every time, snap-on sealed top eliminates contamination from dust or dirt

LxW (mm)	Packed	Pk	Cat. No.
152x228	100 wipes/container	100	115-0128

Wipe rolls for impregnation, with dispensing bucket

These strong centrefeed wipes made from 70% viscose and 30% polyester are supplied with a resealable dispensing bucket for custom impregnation.

- High absorption capacity
- For pre-impregnation at the point of use
- Can be impregnated with any solution



Description	Colour	LxW (mm)	Packed	Pk	Cat. No.
6 centrefeed rolls and 1 dispensing bucket	White	360x270	6 rolls (90 wipes each) + 1 bucket	6 Roll	115-2420

Kitchen wipes, HACCP



These reusable wipes are made from 55% viscose, 35% polyester and 10% binder. Available in four different colours.

- Can be rinsed and reused
- High absorption
- Non scratching
- Recyclable
- Food contact certified
- HACCP colour coding

Typical applications include: Food & beverage, healthcare, janitorial tasks.

Basis weight: 45 g/m²
Absorption: ≥500%

Description	Colour	LxW (mm)	Packed	Pk	Cat. No.
1/4-folded	Blue	380x600	6 packs of 50 wipes	300	115-2400
1/4-folded	Green	380x600	6 packs of 50 wipes	300	115-2402
1/4-folded	Yellow	380x600	6 packs of 50 wipes	300	115-2404
1/4-folded	Red	380x600	6 packs of 50 wipes	300	115-2413

Industrial wipes, clean smooth, turquoise

These spun-lace wipes, made from 55% wood pulp and 45% polyester, are extremely strong, even when wet.

- Low fibre and particle release
- Resistant to acids and solvents
- Good absorption
- Do not contain any binders, silicone-free
- Recyclable
- Food contact certified

Typical applications include: Pharma, aerospace, semiconductor, food & beverage, automotive, laboratories, printing and healthcare.


Basis weight: 73 g/m²

Absorption: 600%



Description	Colour	LxW (mm)	Packed	Pk	Cat. No.
Industrial roll (500 wipes)	Turquoise	380x300	2 rolls, 500 wipes each	2 Roll	115-2396
Z-folded	Turquoise	380x300	10 packs, 50 wipes each	500	115-2397

Industrial wipes, clean smooth, white



These bleached spunlace wipes made from 55% wood pulp and 45% polyester are extremely strong, even when wet.

- Low fibre and particle release
- Resistant to acids, solvents and chlorine
- Good absorption
- Do not contain any binders, silicone-free
- Recyclable
- Food contact certified

Typical applications include: Pharma, aerospace, semiconductor, food & beverage, automotive, laboratories, printing, healthcare.

Basis weight: 68 g/m²

Absorption: 400%

Description	Colour	LxW (mm)	Packed	Pk	Cat. No.
Industrial roll (500 wipes)	White	380x300	2 rolls, 500 wipes each	2 Roll	115-2394
Z-folded	White	380x300	10 packs, 50 wipes each	500	115-2395

Industrial wipes, raspy

These spunlace apertured wipes, made from 51% wood pulp and 49% polyester, are slightly abrasive on one side and soft on the other.

- Extremely strong
- Slightly abrasive side for heavy-duty action
- Resistant to acids, solvents and chlorine
- Recyclable
- Food contact certified

Typical applications include: Food & beverage, automotive, industry, printing, laboratories, healthcare, janitorial tasks.

Basis weight: 75 g/m²

Absorption: 760%



Description	Colour	LxW (mm)	Packed	Pk	Cat. No.
Industrial roll (400 wipes)	White	380x270	1 roll with 400 wipes	1 Roll	115-2398
Z-folded	White	380x270	10 packs, 35 wipes each	350	115-2399

Industrial wipes, melt-blown

These melt-blown wipes, made from 100% PP, have a very high absorption capacity.

- Absorb water, grease and oil
- Resistant to acids and solvents
- Soft, synthetic, colourfast
- Recyclable
- Silicone-free
- Food contact certified

Typical applications include: Aerospace, automotive, electronics, food & beverage, laboratories, printing and janitorial tasks.


Basis weight: 70 g/m²

Absorption: 660% (water); 1350% (oil)



Description	Colour	LxW (mm)	Packed	Pk	Cat. No.
Industrial roll (500 wipes)	Blue	380x320	1 roll with 500 wipes	1 Roll	115-2418
Z-folded	Blue	380x400	20 packs, 25 wipes each	500	115-2419

Industrial wipes, DRC63



These wet-laid, double re-creped, general purpose wipes, made from wood pulp with binder, can be used for many different applications.

- Excellent absorption capacity of water and oil
- Resistant and non scratching
- Stronger and more absorbent than standard kitchen rolls
- Recyclable
- Food contact certified

Typical applications include: Food & beverage, automotive, industry, printing, healthcare and janitorial tasks.

Basis weight: 63 g/m²

Absorption: 760%

Description	Colour	LxW (mm)	Packed	Pk	Cat. No.
Industrial roll (700 wipes)	White	350x260	1 roll with 700 wipes	1 Roll	115-2392
Z-folded	White	380x280	20 packs, 50 wipes each	1.000	115-2393

Industrial wipes, super soft


These very soft wipes, made from a spun-lace fabric, consisting of 60% viscose and 40% polyester, are ideal for sensitive surfaces and light-duty cleaning.

- Non scratching: Ideal for polishing tasks
- Very absorbent
- Resistant to bases and solvents
- Recyclable, does not contain any binders


Typical applications include: Automotive, electronics, optics, industry, laboratories, healthcare and janitorial tasks.

Weight: 50 g/m²

Absorption 1010%



Description	Colour	LxW (mm)	Packed	Pk	Cat. No.
Industrial roll (375 wipes)	White	400x300	2 rolls, 375 wipes each	2 Roll	115-2416
Z-folded	White	380x320	10 packs, 50 wipes each	500	115-2417



Cleanroom wipes, Spec-Wipe® 3E

Durable, non woven wipes, made from 46% polyester and 54% cellulose.

- Hydroentangled: For a durable wipe that produces less particulation than some other non woven methods
- Hydroentangled without chemicals to reduce ionic contamination
- Lightweight polyester/cellulose combination: Economical and sorbent
- Creped fabric for softness and bulk

Suitable for use in Class 1000, ISO 6, M4.5 cleanrooms.

Colour: White

Model	Spec-Wipe® 3E
Absorbency - extrinsic (ml/m ²)	283
Absorbency - intrinsic (ml/g)	5,43
Absorbency - rate (s)	1
Basis weight (g/m ²)	52
Extractables - using deionised water (g/m ²)	0,025
Extractables - using isopropyl alcohol (g/m ²)	0,020
Fibres (>100 µm) (×10 ³ /m ²)	44
Ions - chloride (ppm)	9,4
Ions - sodium (ppm)	8,3
Particle generation (>=0.5 µm) (×10 ⁶ /m ²)	58,2

L×W (mm)	Packed	Pk	Cat. No.
229×229	300 per bag	300	115-0597
305×305	150 per bag	150	115-0598

Dispenser for wipe rolls



These dispensers, made from white lacquered metal, are suitable for use with wipe rolls.

- Floor standing or wall mounting
- For wipe rolls with a max. outer diameter of 40 cm and a max. width of 42 cm

Description	Pk	Cat. No.
Dispenser, white lacquered metal, floor standing, for wipe rolls (max. Ø 400 mm, max. width 420 mm)	1	115-2424
Dispenser, white lacquered metal, wall mounting, for wipe rolls (max. Ø 400 mm, max. width 420 mm)	1	115-2425

Cleaning solution, IPASEPT 70®

Iso-propyl alcohol Ph. Eur. 70% V/V
with purified water Ph. Eur. 30% V/V,
filtered at 0,2 µm

- Filled under GMP conditions in cleanroom ISO Class 7
- Certificate of Analysis is available for each batch
- Fine stream for precise applications, or wide spray for covering large areas

For surface spray cleaning with or without wipes. Also for cleaning gloved hands before entry into cleanrooms. Ideal for wiping down entrances to controlled environments.





Description	Packed	Pk	Cat. No.
IPASEPT 70®	6×1 l trigger spray bottles	6 l	93002.1016

Silicone grease for lubrication



Description	Pk	Cat. No.
Silicone grease TECHNICAL	50 g	6674.0050

Cleanroom anti-fatigue mats, bubble and smooth



These mats offer the comfort of a polyurethane mat and is suitable for cleanroom ISO Class 5 (Class 100) environments. It is either available with a bubble surface, that features a unique ergonomic and anti-fatigue design, or with a smooth surface that is ideal for areas that combine sitting and standing.

- Extremely good chemical and UV resistance (except for organic solvents)
- Free from silicone and latex
- Closed cell / non porous surface allows for easy wash down and cleaning
- Smooth mat: Ideal if carts or trolleys must be used
- Bubble mat: Ergonomic benefit of stimulating the body to increase blood circulation
- Warranty up to 12 years


DIN 53505, EN ISO 14664-1



115-0745

LxW (mm)	Thickness (mm)	Colour	Pk	Cat. No.
Smooth surface, blue				
600x900	11	Blue	1	115-0735
600x1800	11	Blue	1	115-0736
900x1200	11	Blue	1	115-0739
900x1800	11	Blue	1	115-0740
1200x1800	11	Blue	1	115-0742
Smooth surface, grey				
600x900	11	Grey	1	115-0745
600x1800	11	Grey	1	115-0746
900x1200	11	Grey	1	115-0749
900x1800	11	Grey	1	115-0756
1200x1800	11	Grey	1	115-0758
Bubble surface, blue				
600x900	14,5	Blue	1	115-0715
600x1200	14,5	Blue	1	115-1061
600x1800	14,5	Blue	1	115-0716
900x1200	14,5	Blue	1	115-0719
900x1800	14,5	Blue	1	115-0720
1200x1200	14,5	Blue	1	115-1090
1200x1800	14,5	Blue	1	115-0722
Bubble surface, grey				
600x900	14,5	Grey	1	115-0725
600x1200	14,5	Grey	1	115-1106
600x1800	14,5	Grey	1	115-0726
900x1200	14,5	Grey	1	115-0729
900x1800	14,5	Grey	1	115-0730
1200x1200	14,5	Grey	1	115-1135
1200x1800	14,5	Grey	1	115-0732

Ergonomic mats, basic bubble and basic smooth



115-1149



115-1202

These mats are made of low density polyurethane providing an ergonomic, anti-fatigue profile. Bubble mats are designed for dry, non critical environments, reducing trip hazards, body aches, quality issues and sick leave. These smooth, flat mats are ideal for dry, non critical environments where sitting and standing work areas are combined, or where carts or trolleys must be used.

- Moulded bevelled edges will not curl up, thus eliminating trip hazards
- Any size or any shape - will always be one piece: No overlapping or worn out connecting mats
- Ergonomic design stimulates the body to increase blood circulation and alertness
- Silicone- and latex-free
- Warranty up to 3 years

DIN 53505, EN ISO 527-1, ASTM D 624

LxW (mm)	Thickness (mm)	Colour	Pk	Cat. No.
Basic Bubble				
600x900	15	Grey	1	115-1149
600x1200	15	Grey	1	115-1150
600x1500	15	Grey	1	115-1151
600x1800	15	Grey	1	115-1152
600x2100	15	Grey	1	115-1153
600x2400	15	Grey	1	115-1154
600x2700	15	Grey	1	115-1155
600x3000	15	Grey	1	115-1156
600x3300	15	Grey	1	115-1157
600x3600	15	Grey	1	115-1158
600x3900	15	Grey	1	115-1159
600x4200	15	Grey	1	115-1160
600x4500	15	Grey	1	115-1161
600x4800	15	Grey	1	115-1162
600x5100	15	Grey	1	115-1163
600x5400	15	Grey	1	115-1164
600x5700	15	Grey	1	115-1165
600x6000	15	Grey	1	115-1166
900x900	15	Grey	1	115-1167
900x1200	15	Grey	1	115-1168
900x1500	15	Grey	1	115-1169
900x1800	15	Grey	1	115-1170
900x2100	15	Grey	1	115-1171
900x2400	15	Grey	1	115-1172
900x2700	15	Grey	1	115-1173
900x3000	15	Grey	1	115-1174
900x3300	15	Grey	1	115-1175
900x3600	15	Grey	1	115-1176
900x3900	15	Grey	1	115-1177
900x4200	15	Grey	1	115-1178
900x4500	15	Grey	1	115-1179
900x4800	15	Grey	1	115-1180
900x5100	15	Grey	1	115-1181
900x5400	15	Grey	1	115-1182
900x5700	15	Grey	1	115-1183
900x6000	15	Grey	1	115-1184
1200x1200	15	Grey	1	115-1185
1200x1500	15	Grey	1	115-1186

LxW (mm)	Thickness (mm)	Colour	Pk	Cat. No.
Basic Bubble				
1200x1800	15	Grey	1	115-1187
1200x2100	15	Grey	1	115-1188
1200x2400	15	Grey	1	115-1189
1200x2700	15	Grey	1	115-1190
1200x3000	15	Grey	1	115-1191
1200x3300	15	Grey	1	115-1192
1200x3600	15	Grey	1	115-1193
1200x3900	15	Grey	1	115-1194
1200x4200	15	Grey	1	115-1195
1200x4500	15	Grey	1	115-1196
1200x4800	15	Grey	1	115-1197
1200x5100	15	Grey	1	115-1198
1200x5400	15	Grey	1	115-1199
1200x5700	15	Grey	1	115-1200
1200x6000	15	Grey	1	115-1201
Basic Smooth				
600x900	11	Grey	1	115-1202
600x1200	11	Grey	1	115-1203
600x1500	11	Grey	1	115-1204
600x1800	11	Grey	1	115-1205
600x2100	11	Grey	1	115-1206
600x2400	11	Grey	1	115-1207
600x2700	11	Grey	1	115-1208
600x3000	11	Grey	1	115-1209
600x3300	11	Grey	1	115-1210
600x3600	11	Grey	1	115-1211
600x3900	11	Grey	1	115-1212
600x4200	11	Grey	1	115-1213
600x4500	11	Grey	1	115-1214
600x4800	11	Grey	1	115-1215
600x5100	11	Grey	1	115-1216
600x5400	11	Grey	1	115-1217
600x5700	11	Grey	1	115-1218
600x6000	11	Grey	1	115-1219
900x900	11	Grey	1	115-1220
900x1200	11	Grey	1	115-1221
900x1500	11	Grey	1	115-1222
900x1800	11	Grey	1	115-1223
900x2100	11	Grey	1	115-1224
900x2400	11	Grey	1	115-1225
900x2700	11	Grey	1	115-1226
900x3000	11	Grey	1	115-1227
900x3300	11	Grey	1	115-1228
900x3600	11	Grey	1	115-1229
900x3900	11	Grey	1	115-1230
900x4200	11	Grey	1	115-1231
900x4500	11	Grey	1	115-1232
900x4800	11	Grey	1	115-1233
900x5100	11	Grey	1	115-1234
900x5400	11	Grey	1	115-1235
900x5700	11	Grey	1	115-1236
900x6000	11	Grey	1	115-1237
1200x1200	11	Grey	1	115-1238
1200x1500	11	Grey	1	115-1239
1200x1800	11	Grey	1	115-1240
1200x2100	11	Grey	1	115-1241
1200x2400	11	Grey	1	115-1242
1200x2700	11	Grey	1	115-1243
1200x3000	11	Grey	1	115-1244
1200x3300	11	Grey	1	115-1245
1200x3600	11	Grey	1	115-1246
1200x3900	11	Grey	1	115-1247
1200x4200	11	Grey	1	115-1248
1200x1800	11	Grey	1	115-1249
1200x4800	11	Grey	1	115-1254
1200x1800	11	Grey	1	115-1255
1200x5400	11	Grey	1	115-1256
1200x5700	11	Grey	1	115-1257
1200x6000	11	Grey	1	115-1258

Austria

VWR International GmbH
Graumannsgasse 7
1150 Wien
Tel.: +43 1 97 002 0
Email: info.at@vwr.com

Belgium

VWR International bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
3001 Leuven
Tel.: +32 (0) 16 385 011
Email: vwr.be@vwr.com

Czech Republic

VWR International s. r. o.
Veetee Business Park
Pražská 442
CZ - 281 67 Stříbrná Skalice
Tel.: +420 321 570 321
Email: info.cz@vwr.com

Denmark

VWR International A/S
Tobaksvej 21
2860 Søborg
Tel.: +45 43 86 87 88
Email: info.dk@vwr.com

Finland

VWR International Oy
Valimotie 9
00380 Helsinki
Tel.: +358 (0) 9 80 45 51
Email: info.fi@vwr.com

France

VWR International S.A.S.
Le Périgares – Bâtiment B
201, rue Carnot
94126 Fontenay-sous-Bois cedex
Tel.: 0 825 02 30 30* (national)
Tel.: +33 (0) 1 45 14 85 00 (international)
Email: info.fr@vwr.com
* 0,18 € TTC/min

Germany

VWR International GmbH
Hilpertstraße 20a
D - 64295 Darmstadt
Tel.: 0800 702 00 07* (national)
Tel.: +49 (0) 6151 3972 0 (international)
Email: info.de@vwr.com
*Freecall

Hungary

VWR International Kft.
Simon László u. 4.
4034 Debrecen
Tel.: +36 52 521130
Email: info.hu@vwr.com

Ireland / Northern Ireland

VWR International Ltd / VWR International
(Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: +353 (0) 1 88 22 222
Email: sales.ie@vwr.com

Italy

VWR International S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: +39 02 3320311
Email: info.it@vwr.com

The Netherlands

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: +31 (0) 20 4808 400
Email: info.nl@vwr.com

Norway

VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: +47 22 90 00 00
Email: info.no@vwr.com

Poland

VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: +48 58 32 38 200
Email: info.pl@vwr.com

Portugal

VWR International - Material de
Laboratório, Lda
Centro Empresarial de Alfragide
Rua da Indústria, nº 6
2610-088 Alfragide
Tel.: +351 21 3600 770
Email: info.pt@vwr.com

Spain

VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: +34 902 222 897
Email: info.es@vwr.com

Sweden

VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: +46 (0) 8 621 34 00
Email: kundservice.se@vwr.com

Switzerland

VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: +41 (0) 44 745 13 13
Email: info.ch@vwr.com

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard - Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: +44 (0) 800 22 33 44
Email: uksales@vwr.com

China

VWR International China Co., Ltd.
Shanghai Branch
Room 256, No. 3058 Pusan Road
Pudong New District
Shanghai 200123
Tel.: +86 21 5898 6888
Email: info_china@vwr.com

India

VWR Lab Products Private Limited
No.139, BDA Industrial Suburb,
6th Main, Tumkur Road, Peenya Post,
Bangalore, India – 560058
Tel.: +91 80 28078400
Email: vwr_india@vwr.com

Singapore

VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel: +65 6505 0760
Email: sales.sg@vwr.com

GO TO **VWR.COM** FOR THE LATEST
NEWS, SPECIAL OFFERS AND
DETAILS OF YOUR LOCAL VWR
DISTRIBUTOR