





MEET REBEL

THE MICROSCOPE TO CHALLENGE THE STATUS QUO.

- REDEFINED USABILITY
- REDESIGNED ERGONOMICS
- REIMAGINED WORKFLOW







UPRIGHT



INVERTED

PRODUCT SPECIFICATIONS



STAND	CAMERA	DISPLAY	DATA
<p>Hybrid Microscope: Performs both UPRIGHT and INVERTED observation modes</p> <p>Removable display</p> <p>Optional binocular module</p>	<p>Brightfield: 8MP CMOS color camera</p> <p>Repeatable camera and LED light settings</p> <p>Anti-Shake feature: Leverages built in accelerometer to minimize vibrations prior to capture</p>	<p>Apple iPad Pro: 12.9" Retina® display</p> <p>2732 x 2048 (264 pixels/inch)</p> <p>Full circular FOV (field of view) as seen through conventional eyepieces</p> <p>Touchscreen interface</p>	<p>Transfer via: USB</p> <p>AirDrop® (iOS and Mac)</p> <p>Seamless integration with file sharing</p>
ERGONOMICS	SOFTWARE	OPTICS	<p>VWR.COM</p> <p>Prices and product details are current when published and subject to change without notice. Certain products may be limited by country, federal, state, provincial, or local regulations. VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. Offers valid in countries listed above, void where prohibited by law or company policy, while supplies last. Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.</p>
<p>Remote Focus Controller for coarse / fine focus</p> <p>TRUE ERGO TABLE by Labconco: Adjustable height table with flush mounted stage for improved workflow</p> <p>REBEL Logic+ Biosafety Cabinet by Labconco: Flush mounted stage for improved workflow</p>	<p>Touch interface for microscope and camera control</p> <p>Measurement annotations: Length, Area, Cell Count, Scale Bar</p> <p>Touch-based cell counting tool (up to 6 different cell types)</p> <p>Downloads and updates available via App Store</p>	<p>1.25x – 100x Plan Achromat and Plan Apochromat options</p> <p>Brightfield, Darkfield, and Phase Contrast capabilities</p> <p>Extra Long Working Distance (ELWD) condenser available</p> <p>High Resolution (High N.A.) Flip-out condenser available</p>	