


# Avantor's J.T.Baker robotic consumables

Introductory offers!





**TECAN STYLE TIPS**

Robotic tips for Tecan and other OEM workstations

**Type Tecan is compatible with the following robotic instruments:**

- Abbott Diagnostic: m2000
- Adaltis: NEXgen
- AXA Diagnostics: SkyLAB 752
- Hologic: Panther Fusion
- NGNY: AQUA System
- Perkin Elmer (Packard): Multiprobe II, Janus with Multiprobe head
- QIAgen/Corbett: CAS 1200/4200/4800, QIAgility (49002-0106 only)
- Tecan: Genesis, Freedom EVO, MiniPrep, Cavro, Fluent (C-Tray only)



Description	Tip colour	Volume, µL	Rack**	Packaging	Sterile/non-sterile	Case of	Cat. No.	List price £	Promo price £
Filtered	Black	50	S-tray	12x 2 racks, 96 each, blistered	Non-sterile/PCR clean	2304	49006-0106UK	313.00	156.50
Filtered	Black	50	S-tray	12x 2 racks, 96 each, blistered	Sterile*	2304	49006-1106UK	352.00	176.00
Filtered	Black	200	S-tray	12x 2 racks, 96 each, blistered	Non-sterile/PCR clean	2304	49002-0106UK	313.00	156.50
Filtered	Black	200	S-tray	12x 2 racks, 96 each, blistered	Sterile*	2304	49002-1106UK	352.00	176.00
Filtered	Black	1000	S-tray	12x 2 racks, 96 each, blistered	Non-sterile/PCR clean	2304	49003-0106UK	352.00	176.00
Filtered	Black	1000	S-tray	12x 2 racks, 96 each, blistered	Sterile*	2304	49003-1106UK	391.00	195.50
Non-filtered	Black	50	S-tray	12x 2 racks, 96 each, blistered	Non-sterile/PCR clean	2304	49006-0006UK	254.00	127.00
Non-filtered	Black	50	S-tray	18x 10 racks, 96 each, cardboard	Non-sterile	17280	49006-0000UK	1760.00	880.00
Non-filtered	Black	200	S-tray	12x 2 racks, 96 each, blistered	Non-sterile/PCR clean	2304	49002-0006UK	254.00	127.00
Non-Filtered	Black	200	S-tray	18x 10 racks, 96 each, cardboard	Non-sterile	17280	49002-0000UK	1760.00	880.00
Non-filtered	Black	1000	S-tray	12x 2 racks, 96 each, blistered	Non-sterile/PCR clean	2304	49003-0006UK	294.00	147.00
Non-filtered	Black	1000	S-tray	10x 10 racks, 96 each, cardboard	Non-sterile	9600	49003-0000UK	1140.00	570.00

\*Pre-order on Sterile products available.

\*\* C-tray racks also available, please contact us for lead times and pricing.



**HAMILTON STYLE TIPS**

Robotic tips for Hamilton Microlab Star™, Nimbus®, Vantage™ and Microlab Prep



Description	Tip colour	Volume, µL	Packaging	Sterile/non-sterile	Case of	Cat. No.	List price £	Promo price £
Filtered	Black	50	12x 5 racks, 96 each, blistered	Non-sterile/PCR clean	5760	49010-0105UK	479.00	311.35
Filtered	Black	50	12x 5 racks, 96 each, blistered	Sterile*	5760	49010-1105UK	535.00	347.75
Filtered	Black	300	16x 5 racks, 96 each, blistered	Non-sterile/PCR clean	7680	49008-0104UK**	660.00	429.00
Filtered	Black	300	12x 5 racks, 96 each, blistered	Non-sterile/PCR clean	5760	49008-0105UK	496.00	322.40
Filtered	Black	300	12x 5 racks, 96 each, blistered	Sterile*	5760	49008-1105UK	524.00	340.60
Filtered	Black	1000	8x 5 racks, 96 each, blistered	Non-sterile/PCR clean	3840	49009-0105UK	323.00	209.95
Filtered	Black	1000	8x 5 racks, 96 each, blistered	Sterile*	3840	49009-1105UK	380.00	247.00
Non-Filtered	Black	50	12x 5 racks, 96 each, blistered	Non-sterile/PCR clean	5760	49010-0003UK	393.00	255.45
Non-Filtered	Black	300	12x 5 racks, 96 each, blistered	Non-sterile/PCR clean	5760	49008-0003UK	410.00	266.50
Non-Filtered	Black	1000	8x 5 racks, 96 each, blistered	Non-sterile/PCR clean	3840	49009-0005UK	285.00	185.25

Microlab Star™, Nimbus® and Vantage™ are registered trademarks of Hamilton AG.

\*Pre-order on Sterile products available.

\*\*Sale until end of stock.


**J.T.BAKER® SRW MAGNETIC PLATES**

J.T.Baker® Plate SRW (square at top, round at bottom well) magnetic plates are made of high purity PP for efficient sample preparation.

- Standard SBS footprint
- High resistance towards commonly used chemicals, solvents and alcohols
- High stability during centrifugation
- Cylindrical wells enable good sample mixing
- Slightly raised rims reduce the risk of cross-contamination
- Alphanumeric code assists registration and identification of individual samples

J.T.Baker® Plate SRW magnetic system is made for the use with:

- KingFisher™ Flex 96 and Presto; QIAgen: Bio Sprint 96.



Description	Colour	Volume, ml	Sterile/non-sterile	Pack size	Cat. No.	List price £	Promo price £
J.T.Baker® Plate SRW Magnetic	Transparent	0,2	Non-sterile	60	43001-0506UK	181.00	108.60
J.T.Baker® Plate SRW Magnetic	Transparent	2,0	Non-sterile	60	43001-0504UK	240.00	144.00
J.T.Baker® Tip comb	Transparent		Non-sterile	60	43001-0505UK	294.00	176.40

**PLATES**

**J.T.BAKER® PLATE MEDIO**

PP, 96-well storage plate with conical wells.

- Due to the tangency of J.T.Baker® Plate Medio plates, less than 50% stacking space is needed
- Alphanumeric code to assist identification and registration of single samples
- Higher well rims allow an airtight seal
- Conical wells aid sample mixing

J.T.Baker® Plate Medio is perfectly made for limited storage areas (fridge)



Description	Colour	Recommended working volume	Well volume	Pack size	Cat. No.	List price £	Promo price £
J.T.Baker® Medio 96-well plate	Transparent	20 - 500 µl	0,65 ml	40	43001-0102UK	132.00	105.60
J.T.Baker® Medio 96-well plate	Transparent	50 - 1000 µl	1,2 ml	32	43001-0101UK	106.00	84.80
J.T.Baker® Medio 96-well plate	Transparent	50 - 2000 µl	2,3 ml	20	43001-0103UK	75.20	60.16
Capmat for J.T.Baker® Plate Medio	Transparent			20	43001-0104UK	37.60	30.08

**J.T.BAKER® PLATE SQUARE WELL**

These deep square well plates, made from virgin PP, are designed to make the interchanging of plates simpler in automated systems.

They are also ideal for housing both cell and bacterial culture.

- Chemical resistant
- V-bottom well shape maximises liquid removal and aids resuspension
- No inner edges to interfere with magnetic bead collection
- Raised well rims to improve heat sealing
- Conical base aids sample concentration, reconstitution and centrifugation



Description	Colour	Well volume	Pack size	Cat. No.	List price £	Promo price £
J.T.Baker® square 24-well plate	Transparent	10,0 ml	20	43001-0066UK	84.60	67.68
J.T.Baker® square 48-well plate	Transparent	5,0 ml	20	43001-0062UK	84.60	67.68
J.T.Baker® square 96-well plate	Transparent	2,0 ml	20	43001-0020UK	75.20	60.16
Capmat for J.T.Baker® Plate Square Well 24	Transparent		20	43001-0070UK	48.90	39.12
Capmat for J.T.Baker® Plate Square Well 48	Transparent		20	43001-0069UK	43.30	34.64
Capmat for J.T.Baker® Plate Square Well 96	Transparent		20	43001-0022UK	37.60	30.08



# SCIENCE DELIVERED

From scientific discovery to scale-up and commercial delivery, you need mission-critical products, services and solutions on a global scale.

Avantor® offers global access to our own portfolio of trusted, quality brands and critical products through our premier delivery platform, VWR®. All of this, combined with infrastructure strategically located to help serve your specific needs, helps move science forward – fast. That's **science delivered**.

**Get the mission-critical scientific products and solutions you need from Avantor, delivered through VWR.**

➔ [vwr.com/science-delivered](https://vwr.com/science-delivered)



delivered by **VWR**™

---

#### CONTACT:

For more information please contact [uksales@vwr.com](mailto:uksales@vwr.com)